



ViziRev™

Ascending Resettable Temperature Solution

ATI's ViziRev Technology allows users a clear, discernible, visual indication of product temperature status.


Providing an accurate visual indication. ATI's ViziRev Technology is an exceptional way to show temperature change on a variety of products and materials, offering a wide range of flexibility for many industries and markets. The ViziRev Technology can be customized for product specific applications. Transitioning from one image to another during temperature change, ViziRev provides a clear cut indication of temperature fluctuation.



APPLICATIONS INCLUDE:

- Vaccines
- Injectables
- Blood Bags and IV's
- Biologics
- Tissues
- TIC Panels/Controlled Shipping Materials
- Slush Foods
- Specialty chemicals and materials
- In process checks

DIMENSIONS



SPECIFICATIONS

Color change	White to reveal
Temperature capability range	0°C to 18°C
Commonly requested events	5°C, 8°C, 10°C, 18°C
Accuracy	±1°C
Activation specs	Must be chilled a minimum of 10°C below response temperature to activate



RESETTABLE COLOR CHANGING TECHNOLOGY:

- Patent pending
- Primary/Secondary label
- Shelf life data available
- Activatable/Resettable
- Requires no special shipping
- Customizable